

TDS – Prontosan[®] Wound Gel | Prontosan[®] Wound Gel X

Cleansing and moistening of skin wounds and burns. For the prevention of biofilm.



- Removes & prevents biofilm
- Prevents infections
- Reduces healing time
- Painless dressing changes
- Well-known substances with low allergenic potential
- Compatible with commonly used wound dressings
- Can be used up to 8 weeks after first opening

STORAGE CONDITIONS

Store +5 to +30 °C

SHELF-LIFE

36 months

STABILITY

Can be used up to 8 weeks after first opening

While necrotic tissue and slough can be optimally removed by Prontosan[®] Wound Irrigation Solution, it is necessary to first remove heavy encrustations by keeping the wound moist for a longer period of time. For this task, Prontosan[®] Wound Gel and Gel X provide the optimal solution.

WHEN TO USE WHICH GEL

Prontosan [®] Wound Gel	Prontosan [®] Wound Gel X		
30 ml pod	250 g tube		
Fluid	Highly viscous		
For smaller, deep and tunneling wounds, as well as for difficult to access wounds For burn wounds I. & II. degree	For large surface area wounds, e.g. • Ulcus cruris • Decubitus For all burn wounds, incl. III. & IV. degree		
If little quantity is required	When large quantities are required Where Prontosan® Wound Gel is too fluid		

B BRAUN SHARING EXPERTISE

TDS – Prontosan[®] Wound Gel | Prontosan[®] Wound Gel X

Cleansing and moistening of skin wounds and burns. For the prevention of biofilm.

HINTS AND TIPS

All wounds should in principle first be rinsed and cleansed with Prontosan[®] Wound Irrigation Solution, Prontosan[®] Gel and Gel X remain on the wound until the next dressing change. It therefore has a long acting effect.


INDICATIONS


Cleansing, moistening and decontamination of skin wounds and burns

- Traumatic wounds
- Postoperative wounds
- Chronic skin ulcers
 (e.g. venous, diabetic or pressure ulcers)
- Thermal wounds
- Chemical wounds
- (acid- and alkaliinduced)
- Radiation-induced wounds

Moistening of bandages and wound dressings







APPLICATION

The choice between Prontosan[®] Wound Gel and Gel X allows an optimal application in large surface area as well as deep wounds.

- For large surface area wounds apply a 3-4 mm thick layer of the **Prontosan**[®] **Wound Gel X** and cover with a secondary dressing.
- 2 For the application in deep or tunneling wounds, wound cavities and difficult to access areas, apply a 3-5 mm layer of **Prontosan**[®] **Wound Gel** and cover with a secondary dressing.

Prontosan® Wound Gel



	Size	Pcs/Pack	Article number
Prontosan [®] Wound Gel	30 ml	1	400515
Prontosan [®] Wound Gel X	50 g	1	400517
Prontosan [®] Wound Gel X	250 g	1	400508

Composition: 0.1 % Polyaminopropyl Biguanide (Polihexanide), Betaine Surfactant, Purified Water, Glycerol, Hydroxyethylcellulose