



GENERAL
OPEN SURGERY

CSC-KOL[®] INTRALUMINAL STAPLER FOR SINGLE USE

INNOVATION IN LOW COLORECTAL ANASTOMOSIS
WITH TRANSANAL CONTROL

CSC-KOL[®] INTRALUMINAL STAPLER FOR SINGLE USE

INNOVATION IN LOW COLORECTAL ANASTOMOSIS WITH TRANSANAL CONTROL

EXISTING CLINICAL PROBLEMS




Fig. 1: This picture shows the crossings of the staple lines as well as the formation of the dog ears.

- Total mesorectal excision (TME) is the gold standard technique for rectal surgery as it decreases local recurrence and improves functional results¹.
- On the other hand, the introduction of the double stapling technique (DST) has greatly facilitated the anastomosis after low anterior resection (LAR); it obviates the use of lower purse-string suture and permits a lower and easier anastomosis².
- However, stapling across staple lines and the lateral intersecting margins created (so-called dog ears), have been reported to increase the risk of stenosis¹¹ and anastomotic leakages^{3,4}, which is the most common complication after rectal cancer surgery and its rate remains at approximately 10%⁵⁻¹⁰.




Fig. 2




Fig. 3

**KOL® – TECHNIQUE
FOR A REAL TRANSANAL END-TO-END ANASTOMOSIS**

1. Transanally suture the distal staple line

- Firstly stitch the left end of the staple line. Suture another three or four points until the right end of the staple line. (Fig. 2)
- Pull the two groups of sutures into the right and left traction holes separately, by using the suture threader. (Fig. 3)

2. Anastomosis

- Close the stapler into the green zone, while pulling the staple line into the housing and fire the stapler.
- Check the specimen to make sure the complete resection of the distal staple line.

The picture shows the benefits of a real transanal anastomosis without intersecting margins (dog ears) and crossings of staple lines, reducing the risk of leakages and stenosis.





Fig. 4

PRODUCT FEATURES AND BENEFITS

SINGLE USE
INTRALUMINAL
STAPLER
CSC-KOL®
INNOVATION



Power-Cut

Improved technology to cut through intersecting staples.

Double specimen housing

12.6 - 17.8 cm³ to accommodate the staple line and the desired amount of tissue to be resected.

Just-Fit & Uni-Link

enhance direct firing force transition and use parallel tissue compression and staple formation.

Four Traction Holes

allow to pull the staple line and the tissue into the housing.


REFERENCES

- ¹ Kawada K, Hasegawa S, Koya H, Hirai K, Okoshi K, Nomura A, et al. Risk factors for anastomotic leakage after laparoscopic low anterior resection with DST anastomosis. *Surg Endosc.* 2014;28:2988-95.
- ² Cohen Z, Myers E, Langer B, Taylor R, Railton RH, Jamieson C. Double stapling technique for low anterior resection. *Dis Colon Rectum.* 1983;26(4):231-5.
- ³ Kuroyanagi H, Oya M, Ueno M, Fujimoto Y, Yamaguchi T, Muto T. Standardized technique of laparoscopic intracorporeal rectal transection and anastomosis for low anterior resection. *Surg. Endosc.* 2008;22:557-61.
- ⁴ Rudi M, Roumen H, Rahusen F, Wijnen M, Van Uchelen F. "Dog ear" formation after double-stapled low anterior resection as a risk factor for anastomotic disruption. *Dis Colon Rectum.* 2000;43(4):522-5.
- ⁵ Kang CY, Halabi WJ, Chaudhry OO, Nguyen V, Pigazzi A, Carmichael JC, et al. Risk factors for anastomotic leakage after anterior resection for rectal cancer. *JAMA Surg.* 2013;148:65-71.
- ⁶ den Dulk M, Marijnen C, Collette L, Putter H, Pahlman L, Folkesson J, et al. Multicentre analysis of oncological and survival outcomes following anastomotic leakage after rectal cancer surgery. *Br J Surg.* 2009;96:1066-75.
- ⁷ Peeters KC, Tollenaar R, Marijnen C, Kranenbarg E, Steup W, Wiggers T. Risk factors for anastomotic failure after total mesorectal excision of rectal cancer. *Br J Surg.* 2005;92:211-6.
- ⁸ Paun BC, Cassie S, MacLean A, Dixon E, Buie W. Postoperative complications following surgery for rectal cancer. *Ann Surg.* 2010;251:807-18.
- ⁹ Snijders HS, Wouters M, Van Leersum NJ, Kolfsoorten NE, Henneman D, De Vries AC, et al. Meta-analysis of the risk for anastomotic leakage, the postoperative mortality caused by leakage in relation to the overall postoperative mortality. *Eur J Surg Oncol.* 2012;38:1013-9.
- ¹⁰ Park JS, Choi G, Kim S, Kim H, Hyeong R, Kim N, et al. Multi-center analysis of risk factors for anastomotic leakage after laparoscopic rectal cancer excision: the Korean laparoscopic colorectal surgery study group. *Ann Surg.* 2013;257:665-71.
- ¹¹ Griffen F, Knight C, Whitaker J, Knight C. The double stapling technique for low anterior resection. Results, modifications and observations. *Ann Surg.* 1990;211(6):745-52.

ORDER INFORMATION

Article No.	Anvil Diameter	Colour Code	Staple Quantity	Blade Diameter	Staple Height	Closed Staple Height	Housing Length / Volume
CSC29-KOL	29 mm	Blue	24	20.5 mm	5.0 mm	1.0 ~ 2.5 mm	4 cm / 12.6 cm ³
CSC33-KOL	33 mm	Green	32	24.8 mm	5.0 mm	1.0 ~ 2.5 mm	4 cm / 17.8 cm ³

Transanal Accessories



Manufacturer acc. to MDD 93/42/EEC

Touchstone International Medical Science Co., Ltd. | 21 A Science Plaza
International Science Park | No. 1355 Jinjihu Avenue | Suzhou | China

Distributed by:

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain
Phone +34 93 586 6200 | Fax +34 93 699 6330 | www.bbraun.com

AESCULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The product trademark "KOL" is a registered trademark of Touchstone International Medical Science Co., Ltd. Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.