Distributed by / Vertrieb / Distribuido por / Distribué par / Distribuito da / Gedistribueerd door / توزيع
من **لبيق ن** / **Διανέμεται από** / **販売元** / Diedarkan oleh / Distribuert av / Dystrybutor / Distribuido por /
Distributor / Distribuerad av / Distribütör / Izplatītājs / **Дистрибьютор** / Turustaja:

Effective

GB B. Braun Medical Ltd.
Sheffield, S35 2PW - United Kingdom

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












B | BRAUN

1-Piece Range

EN Instructions for use 1–Piece Range**What Proxima®, Proxima®+, Proxima® Border and Almarys® Optima is and what is used for**

Proxima®, Proxima®+, Proxima® Border and Almarys® Optima are 1-piece collection pouches for enterostomates or urostomates either:

- Closed for colostomates with an integrated filter for deodorization and evacuation of gases.
- Drainable for enterostomates with an adhesive soft clamp or reusable rigid clamp and a large evacuation sleeve allowing the evacuation of thick effluents. Drainable pouches are available with or without filter.
- Drainable for urostomates with a drainage tap (H-I-J) allowing connection to a night bag.*

Proxima®, Proxima®+, Proxima® Border and Almarys® Optima may be used concomitantly with B. Braun stoma care accessories.

Intended purpose and indication:

The product is indicated for collection and containment of stools or urine by patients with enterostomy or urostomy. The product is not intended to be used on a skin or mucus membrane presenting a pathological change or change following disease or a wound.

Closed pouches are primarily intended for colostomates; drainable pouches for colostomates with soft or liquid stools and ileostomates; uro pouches for urostomates.

Clinical benefits & performances: Ostomy appliances compensate for the loss of faecal and/or urinary incontinence.

Patient population: enterostomates and urostomates adults of children above four years.

Contraindications:

There are no formal contraindications for the use of Proxima®, Proxima®+, Proxima® Border and Almarys® Optima.

For the closed pouch Proxima®, Proxima®+, we do not recommend using them by colostomates with liquid stools.

Residual risks or side effects:

Peristomal skin complications can occur while using Proxima®, Proxima®+, Proxima® Border and Almarys® Optima.

Leakages through the filter can occur while using drainable pouches.

1. Preparing the skin

Wash your hands before starting the procedure.

Gently clean the peristomal area with lukewarm water. Avoid ether, alcohol or any other topical product which could interfere with skin protector adhesion. Carefully dab dry the skin with soft and non-shedding cloths.

2. Applying the pouch (closed and drainable pouches)

For drainable pouches, close the closure outlet:

- Drainable pouches with adhesive clamp: start by removing the release paper from the adhesive strip (be careful not to place your fingers on the adhesive part). Centre and stick the adhesive clamp at the bottom of the opening (on the outer side of the pouch, not on the skin side). The width between the clamp and the bottom of the outlet should be equivalent to the width of the clamp (E). Fold the extremity of the sleeve four times and the outlet backwards using the clamp (F). Firmly squeeze the extremities of the clamp to ensure complete closure (G).
- Drainable pouches with rigid clamp (for Proxima+ references only): Open the clamp by pressing on the latch. Place the thin curved part of the clamp towards the plastic part of the outlet (K), fold the film upwards then close the clamp (L-M). If you prefer, use a rigid clamp (to be ordered separately: ref. F00727).
- Drainable "uro" pouches: Make sure the outlet is closed by turning clockwise the tap (I).

Cut the central opening of the hydrocolloid adhesive to the size of your stoma (A) possibly with the help of the stoma guide.

MD

- EN Medical device as per 2017-745/EU
- DE Medizinprodukt gemäß 2017-745/EU
- ES Producto sanitario según el Reglamento (UE) 2017-745
- FR Dispositif médical selon 2017-745/EU
- IT Dispositivo medico secondo 2017-745/UE
- NL Medisch hulpmiddel conform 2017-745/EU
- AR جهاز طبي وفقاً للاتحاد الأوروبي 2017-745
- EL Ιατροτεχνολογικό προϊόν σύμφωνα με τον Κανονισμό 2017/745/ΕΕ
- ET Medicinsk udstyr i henhold til 2017-745/EU
- JA 2017-745/EUに従った医療機器
- LV Medicīniskā ierīce saskaņā ar 2017-745/ES
- MS Peranti perubatan berdasarkan 2017-745/EU
- NO Medisinsk utstyr iht. 2017-745/EU
- PL Wyrób medyczny zgodnie z rozporządzeniem (UE) 2017/745
- PT Dispositivo médico de acordo com a diretiva 2017-745/EU
- RU Медицинское изделие согласно 2017-745/EU
- RO Dispozitiv medical conform 2017-745/EU
- SV Medicinteknisk produkt enligt 2017-745/EU
- TR 2017 745/EU uyarınca tibbi cihaz

EN Instructions for use 1–Piece Range

Do not cut beyond the maximum possible diameter and make sure that there is always a 2 to 3 mm margin around your stoma in order not to strangle or hurt the stomal mucosa.

Separate the pouch walls by crumpling it slightly to introduce air into it (B).

Grasp the pouch by the sides (C).

Peel off the release from the skin protector and position the bottom of the flange opening around the base of the stoma (C). Press the bottom part of the skin protector smoothly onto the skin, then place and press onto the upper part taking care not to make folds.

In order to improve the adhesion, place hands over the skin protector to warm it up (D).

3. Emptying of drainable pouches

Position yourself appropriately.

- Drainable pouches with adhesive clamp: slightly lift upwards the bottom of the pouch, unfold the clamp's extremities and unfold the outlet to empty the pouch. Once the emptying is complete, expel, as much as possible, any matter still in the sleeve by pressing it lightly from top to bottom. Carefully clean its extremity avoiding the use of shedding cloths.
- Drainable pouch with rigid clamp: slightly lifts upwards the bottom of the pouch and press the latch (N) to open the clamp. Remove the clamp before emptying the pouch. Then direct the open end of the pouch over the toilet. Once the emptying is complete, expel, as much as possible, any matter still in the sleeve by pressing it lightly from top to bottom. Carefully clean its extremity avoiding the use of shedding cloths.
- Drainable "Uro" pouches: point the tap towards the toilet bowl, slowly turn the tap counterclockwise until the desired flow rate is obtained (H).

After emptying of drainable pouches, close the closure outlet again as indicated above.

4. Using uro connectors

The connection of the pouch to a B. Braun distal bag can require the use of the supplied uro connector. To clip on the uro connector, press the two wings and insert it into the pouch tap (J). Open the tap turning it half a turn counterclockwise to allow urine to flow in the B. Braun distal bag (H). To un-clip, close the tap of the pouch by turning it half clockwise and pull out the uro connector by pressing the two wings.

5. Using filter protection stickers

Before bathing or showering, stick the filter protection sticker over the gas evacuation openings (O) in order to maintain filter performances.

For drainable pouches, in the case of horizontal use of the pouches (night, etc.) stick the filter protection sticker over the gas evacuation openings (O) in order to avoid leakages through the filter. Replace the pouch when a vertical position is recovered.

6. Removing of the pouch

Pull slightly on the removal tab and remove the pouch from the top to the bottom, whilst maintaining gentle pressure on the skin around the stoma with your other hand.

**Availability may vary from one country to another*

If you note any change on your stoma and/or peristomal skin, please contact your health care professional.

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EN **Warning:** Reuse of single-use devices creates a potential risk to the patient. It may lead to contamination and/or impairment of functional capability. Please seek advice from your health care professional before using convex devices. Only the rigid clamp can be re-used (the soft clamp is a single use device). It is advised to always carefully inspect the rigid clamp before using it. The clamp can be used for 1 month, however if you notice a deterioration please use a new one. Ostomy appliances cannot be reprocessed.

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