Distributed by / Vertrieb / Distribuido por / Distribué par / Distribuito da / Gedistribueerd door /  
توزع من قبل / 販売元 / Diedarkan oleh / Distribuert av / Dystrybutor / Distribuido por / **Дистрибутор** /  
Distributor / Distribuerad av / Distribütör / **Διανέμεται από** / Turustaja / Didistribusikan oleh / Izplatītājs /  
**Phân phối bởi** / 经销商

**GB** B Braun Medical Ltd.  
Sheffield S35 2PW - United Kingdom

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**CN** B. Braun Medical (Shanghai) Int'l Trading Co., Ltd.  
Part S,P and Q, No.285 Gang Ao Road, China  
(Shanghai) Pilot Free Trade Zone, Shanghai, China  
Code postal: 200131 Tél: 021-22163000

**RS** Beohem-3 doo  
11000 Belgrade - Serbia

B. BRAUN MEDICAL SAS | 26 Rue Armengaud | 92210 Saint-Cloud | FRANCE

# 2-Piece Range

Effective

Proxima® 2, Proxima® 2+

---


---

Effective

**B | BRAUN**

2-Piece Range

- EN Medical device as per 2017-745/EU
- DE Medizinprodukt gemäß 2017-745/EU
- ES Producto sanitario según el Reglamento (UE) 2017-745
- FR Dispositif médical selon 2017-745/EU
- IT Dispositivo medico secondo 2017-745/UE
- NL Medisch hulpmiddel conform 2017-745/EU
- AR 2017-745 جهاز طبي وفقًا للاتحاد الأوروبي
- EL Ιατροτεχνολογικό προϊόν σύμφωνα με τον Κανονισμό 2017/745/ΕΕ
- ET Medicinsk udstyr i henhold til 2017-745/EU
- ID Läkinnällinlinen laite direktiivin 2017-745/EU mukaisesti
- JA 2017-745/EUに従った医療機器
- LV Mediciniskā ierīce saskaņā ar 2017-745/ES
- MS Peranti perubatan berdasarkan 2017-745/EU
- NO Medisinsk utstyr iht. 2017-745/EU
- PL Wyrób medyczny zgodnie z rozporządzeniem (UE) 2017/745
- PT Dispositivo médico de acordo com a diretiva 2017-745/EU
- RU Медицинское изделие согласно 2017-745/EU
- RO Dispozitiv medical conform 2017-745/EU
- SV Medicinteknisk produkt enligt 2017-745/EU
- TR 2017 745/EU uyarınca tıbbi cihaz
- VI Thiết bị y tế theo quy định 2017-745/EU
- ZH 符合 2017-745/EU 的医疗器械

What are Proxima® 2, Proxima® 2+ and what are used for:

Proxima® 2, Proxima® 2+ are 2-piece collection systems for enterostomates or urostomates comprising of:

- A flat base plate with:
  - a coupling ring for attaching the pouch,
  - a locking system,
  - loops to attach B. Braun Belt,
  - a hydrocolloid adhesive;
- A convex base plate with:
  - a coupling ring for attaching the pouch
  - a locking system
  - loops to attach B. Braun Belt
  - a convex ring
  - a hydrocolloid adhesive

**Different type of pouches:**

- Closed for colostomates with an integrated filter for deodorization and evacuation of gases,
- Drainable for enterostomates with an adhesive soft clamp or reusable rigid clamp and a large evacuation sleeve allowing the evacuation of thick effluents. Drainable pouches are available with or without filter. Adhesive clamps and a rigid clamp (for Proxima® 2+ references only) are provided in each box,
- Drainable for urostomates with a drainage tap and uro connector (L-M-N) \* allowing connection of a night bag.

Proxima® 2, Proxima® 2+ may be used concomitantly with B. Braun stoma care accessories.

Intended purpose and indication:

The product is indicated for collection and containment of stools or urine by patients with enterostomy or urostomy. The product is not intended to be used on a skin or mucus membrane presenting a pathological change or change following disease or a wound.

Closed pouches are primarily intended for colostomates; drainable pouches for colostomates with soft or liquid stools and ileostomates; uro pouches for urostomates.

Clinical benefits & performances: Ostomy appliances compensate for the loss of faecal and/or urinary incontinence.

Patient population: Enterostomates and urostomates adults or

children above four years.

Contraindications:

For the closed pouch Proxima 2® and Proxima® 2+, we do not recommend using them by colostomates with liquid stools.

Residual risks or side effects:

Peristomal skin complications can occur while using Proxima® 2 and Proxima® 2+.

Leakages through the filter can occur while using drainable pouches.

**1. Preparing the skin**

- Wash your hands before starting the procedure.
- Gently clean the peristomal area with lukewarm water. Avoid ether, alcohol or any other topical product which could interfere with base plate adhesion. Carefully dab dry the skin with soft and non-shedding cloths.

**2. Preparing the base plate**

If you are using a Cut-To-Fit product, cut the central opening of the hydrocolloid adhesive to the size of your stoma (A) possibly with the help of the stoma guide. Do not cut beyond the maximum possible diameter and make sure that there is always a 2 to 3 mm margin around your stoma in order not to strangle or hurt the stomal mucosa.

**3. Applying the base plate**

- Peel off the release.
- Apply the base plate over your stoma in a way such that the locking system is located at the top. In order to improve the adhesion, place hands over the base plate to warm it up (B).
- Ensure the locking system is opened by moving the lever gently upwards without forcing (C).

**4. Applying the pouch: (closed and drainable pouches)**

For drainable pouches, close the closure outlet:

- Drainable pouches with adhesive clamp: start by removing the release paper from the adhesive strip. (be careful not to place your fingers on the adhesive part). Centre and stick the adhesive clamp at the bottom of the opening (on the outer side of the pouch, not on the skin side). The width between the clamp and the bottom of the outlet should be equivalent to the width of the clamp (I). Fold the extremity of the sleeve four times and the outlet backwards using the clamp (J). Firmly

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squeeze the extremities of the clamp to ensure complete closure (K).

- Drainable pouches with rigid clamp (for Proxima 2+ references only) : open the clamp by pressing on the latch. Place the thin curved part of the clamp towards the plastic part of the outlet (O), fold the film upwards then close the clamp (P-Q).
- Drainable "uro" pouches : make sure the outlet is closed by turning clockwise the tap (L).

Separate the pouch walls by crumpling it slightly to introduce air into it.

Grasp the pouch by the sides.

Press the pouch and base plate together starting at the bottom.

Run your fingers around the coupling on the outside of the pouch from the bottom to the top: when a slight click is felt, the pouch is securely positioned on the base plate (D). Push the lever down to lock the system (E). Check that the pouch is securely locked onto the base plate by gently pulling on the pouch.

#### 5. Emptying of drainable pouches

Position yourself appropriately.

- Drainable pouches with adhesive clamp : slightly lift upwards the bottom of the pouch, unfold the clamp's extremities and unfold the outlet to empty the pouch. Once the emptying is complete, expel, as much as possible, any matter still in the sleeve by pressing it lightly from top to bottom. Carefully clean its extremity avoiding the use of shedding cloths.
- Drainable pouch with rigid clamp : slightly lift upwards the bottom of the pouch and press the latch (R) to open the clamp. Remove the clamp before emptying the pouch. Then direct the open end of the pouch over the toilet. Once the emptying is complete, expel, as much as possible, any matter still in the sleeve by pressing it lightly from top to bottom. Carefully clean its extremity avoiding the use of shedding cloths.
- Drainable "Uro" pouches : point the tap towards the toilet bowl, slowly turn the tap counter-clockwise until the desired flow rate is obtained (M).

After emptying of drainable pouches, close the closure outlet again as indicated above.

#### 6. Removing the product

- Removing the pouch : make sure that the tab on the bag is

aligned with the lever on the base plate (F). Open the locking mechanism by moving the locking lever gently upwards (G). Hold the tab on the plastic ring of the pouch then gently pull the pouch out of the base plate, starting at the top and working downwards (H), while maintaining the base plate with your other hand.

- Removing the base plate: pull slightly on the top of the base plate and remove the base plate from the top to the bottom, while maintaining gentle pressure on the skin around the stoma with your other hand.

#### 7. Using uro connectors

The connection of the pouch to a B. Braun distal bag can require the use of the supplied uro connector. To clip on the uro connector, press the two wings and insert it into the pouch tap (N). Open the tap turning it half a turn counterclockwise to allow urine to flow in the B. Braun distal bag. To un-clip, close the tap of the pouch by turning it half clockwise and pull out the uro connector by pressing the two wings (N).

#### 8. Using B. Braun belt

If you wish to use a B. Braun belt, use the two loops on each side.

#### 9. Using filter protection stickers

Before bathing or showering, stick the filter protection sticker over the gas evacuation openings (S) in order to maintain filter performances.

For drainable pouches, in the case of horizontal use of the pouches (night, etc.) stick the filter protection sticker over the gas evacuation openings (S) in order to avoid leakages through the filter. Replace the pouch when a vertical position is recovered.

If you note any change on your stoma and/or peristomal skin, please contact your health care professional.

National regulations for the destruction of potentially infectious materials must be observed.

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**EN Warning:** Reuse of single-use devices creates a potential risk to the patient. It may lead to contamination and/or impairment of functional capability. Please seek advice from your health care professional before using convex devices.

Only the rigid clamp can be re-used (the soft clamp is a single use device). It is advised to always carefully inspect the rigid clamp before using it. The clamp can be used for 1 month, however if you notice a deterioration please use a new one. Ostomy appliances cannot be reprocessed.

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