

Solofix®

Sterile single use steel lancets for capillary blood collection



Solofix®		
REF	200 pieces	6182003
	500 pieces	6185002
Size	Product dimension	39.1 mm x 5.6 mm
	Blade	2.9 mm exposed length
		12° asymmetric point, 18° angle
Material		Stainless steel 1.4016 (X6Cr17) acc. to DIN EN 10088-1
		Compliant with the release limits for nickel as stated in the REACH guidelines (EC) 1907/2006, annex XVII, entry 27
Logistics	Packaging	Single lancet in sterile peel package, 5 pieces per slip, dispenser box with 200/500 pieces
	Shelf life	5 years
Sterilization		Gamma ray irradiation
Certificates		EN ISO 13485
Product classification		Medical Device according to MDD 93/42/EEC, Class IIa

All data refer to typical single values and may be subject to alterations.