

Urimed® Cath

Indwelling Foley Catheters



Wearing time: up to 90 days

Indication: bladder drainage of the urine

- 100 % silicone
- Latex free
- DEHP free


Different lengths available:

- Urimed® Cath Foley Nelaton is 31 cm* (for CH8 and CH10) and 40 cm* (for other CH sizes)
- Urimed® Cath Foley Tiemann is 42 cm*

Different packaging available:


- 5 units boxes with individual packaging
- 10 units boxes

Radiopaque tip and shaft



Silicone catheter material

- Silicone catheters are very gentle for the tissue and hypoallergenic.¹
- Silicone catheters might be preferable to other catheter materials to reduce the risk of encrustation in long-term catheterized patients.¹
- Silicone catheters have a relatively large lumen and have a reduced tendency to encrustation because they are uncoated.¹



* Total length

Sources:

¹ Geng V et al. Evidence-based Guidelines for Best Practices in Urological Health Care: Catheterisation - Indwelling catheters in adults - Urethral and suprapubic. EAUN, 2012

Easy-to-open packaging


- Sterile
- Double packaging
- Two ways to open: lateral or longitudinal
- No-touch design



Tear off the primary peel-pack with the side pre-cuts



and fold it back



or use the pre-cuts to tear it off down the catheter.

Product	Useful length	CH	Ø	Balloon capacity (ml)	10 units		Individual packaging	
					Units/box	Ref.	Units/box	Ref.
 Urimed® Cath Foley Nelaton	23 cm	08	2.7	3	10	4563308	5	4563308-01
		10	3.3	3	10	4563310	5	4563310-01
	33 cm	12	4.0	10	10	4563312	5	4563312-01
		14	4.7	10	10	4563314	5	4563314-01
		16	5.3	10	10	4563316	5	4563316-01
		18	6.0	10	10	4563318	5	4563318-01
		20	6.7	10	10	4563320	5	4563320-01
		22	7.3	10	10	4563322	5	4563322-01
 Urimed® Cath Foley Tiemann	35 cm	12	4.0	10	10	4563412	5	4563412-01
		14	4.7	10	10	4563414	5	4563414-01
	16	5.3	10	10	4563416	5	4563416-01	
	18	6.0	10	10	4563418	5	4563418-01	
	20	6.7	10	10	4563420	5	4563420-01	

B. Braun Avitum AG | 34212 Melsungen | Germany
 Phone +49 5661 71-0 | www.bbraun.com

This document, its contents, including institutional data, information, trademarks and logos mentioned herein are the exclusive property of B. Braun. Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun, is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun. Non binding documents and photographs. For healthcare professional use only.