

Softa® Swabs | Technical Data Sheet

Swabs for Skin Cleansing – Isopropyl Alcohol 70% (non sterile)



	Softa [®] Swabs			
REF	19579	19579		
SWAB SIZE	Folded	32.5 mm x 30 mm		
	Unfolded	65 mm x 30 mm		
INGREDIENTS (INCI NAME) Isopropyl Alcohol, Aqua	ohol, Aqua		
MATERIALS/COMPOSITION	Solution of isopropanol	2-Propanol	70 percent by volume	
		Purified water	ad. 100 volume percent	
		Each swab is soaked with 0.45 g ($ riangle$ 0.5 ml) of an aqueous solution of isopropanol (70 vol. %) equivalent to 0.28 g ($ riangle$ 0.35 ml) isopropanol		
	Wipe	Polypropylene	100 percent by weight	
		No chemical binders		
LOGISTICS	Packaging	Packaging Individually foil-wrapped, air-tight seal Box with 100 pieces		
	Shelf life	4 years		
TERILIZATION	Non sterile			
APPLIED STANDARDS	Regulation EC 1223/200	Regulation EC 1223/2009 on cosmetic products		
	ISO 22716:2007 Cosmet	ISO 22716:2007 Cosmetic guidelines on good manufacturing practices (GMP)		
	Council Directive EC No.	Council Directive EC No. 1907/2006 (SDS)		
	ISO 9001			
PRODUCT CLASSIFICATIO	Cosmetic product			