



GENERAL OPEN SURGERY

AUTO PURSE STRING CLAMP

SINGLE USE INTERNAL STAPLER

AUTO PURSE STRING CLAMP


SINGLE USE INTERNAL STAPLER



The Auto Purse String clamp is a single-use instrument for the placement of temporary purse-string closures in gastric, esophageal, intestinal and colorectal surgery.

It places a circumferential strand of 2–0 nylon suture material, held in place by stainless steel staples.

The stainless steel staples are positioned in the tissue where the purse-string is needed. The device is activated by squeezing the ring handles.



Conventional Staple

HUG Staple

ARTICLE NO.	HANDLE LENGTH	CARTRIDGE LENGTH	SUTURE LENGTH	STAPLE QUANTITY
APS45	310 mm	47 mm	690 mm	18
APS65	256 mm	67 mm	730 mm	28

Distributed by:

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain Phone +34 93 586 6200 | Fax +34 93 588 1096 | www.bbraun.com

AESCULAP® - a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

Manufacturer acc. to MDD 93/42/EEC

Touchstone International Medical Science Co., Ltd. | 21 A Science Plaza International Science Park | No. 1355 Jinjihu Avenue | Suzhou | China

The main product trademark "Aesculap" is a registered trademark of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.