

B. Braun Melsungen AG
 Division Hospital Care
 Marketing & Sales
 34209 Melsungen, Germany

March, 2022

List of Supported Materials KIRO® Oncology

Dear Sir or Madam,

We hereby confirm that the following Products, belonging to the B. Braun Melsungen AG portfolio are currently supported by Kiro Grifols' Automated Drug Admixture solution KIRO® Oncology.

Product Category	Supplier	Product Name	Comments
IV Containers	B. Braun	Ecoflac® plus (50ml, 100ml, 250ml, 500 ml, 1.000 ml)*	Requires the use of P/N 100889 support when used in combination with Secondary Administration Set.
Syringes (Luer Lock)	B. Braun	Omnifix®-F Luer Lock Solo 1 mL 9167006 V Omnifix® Luer Lock Solo 5 mL 4617053 V 10 mL 4617100 V 20 mL 4617207 V 50 mL 4617509 F	Can be used as transfer syringes (except 1ml) or as final containers.
Needles	B. Braun	Sterican® 20 G x 1" REF 4657500	-
Secondary Administration Sets	B. Braun	ProSet Cyto-Set® Mix with needle-based Injection Port and vented spike - REF 4182706	Can be connected to Ecoflac® plus IV container and primed with diluent solution before automated compounding.
Secondary Administration Sets	B. Braun	Cyto-Set® Mix with needle-free access port and vented spike - REF A2907N	Can be connected to Ecoflac® plus IV container and primed with diluent solution before automated compounding. To be used with additional adaptor for automated needle-based injection

Product Category	Supplier	Product Name	Comments
Secondary Administration Sets	B. Braun	Cyto-Set® Mix with needle-free access port, vented spike and 0.2 µm filter - REF A2903N	Can be connected to Ecoflac® plus IV container and primed with diluent solution before automated compounding. To be used with additional adaptor for automated needle-based injection
Secondary Administration Sets	B. Braun	Cyto-Set® Mix with needle-free access port, vented spike and light protected tubing - REF A2908N	Can be connected to IV containers and primed with diluent solution before automated compounding. To be used with additional adaptor for automated needle-based injection
Elastomeric pumps	B. Braun	Easypump® II ST/LT, all references	To be used with additional adaptor for needle-based automated drug admixture in the final container.
Dispensing Pins	B. Braun	Mini-Spike® Chemo V REF 4550587	Used for the reconstitution of powder drug vials.
Syringe Caps	B. Braun	Tamper Evident Syringe Caps REF 418004	Used to deliver automatically capped and ready-to-administer syringes.
Syringe Caps	B. Braun	Puresite Closed Male Connector – REF 4551105 Puresite Dust Cover cap – REF 4551101	Used to deliver automatically capped and ready-to-administer syringes. Non-removable, can be connected to female luer access port.
Irrigation bags	B. Braun	Sterile Water for Irrigation in Ecobag® CE, 3.000 ml REF FREU912	Used for the automated self-cleaning process of KIRO® Oncology.
EVA Mixing bags	B. Braun	Nutrimix® (200 ml, 500 ml, 1.000 ml) REF 2112145, 2112146, 2112147	Used only as empty containers for the collection of discarded diluent.

We hope this information is helpful and comprehensive.

Do not hesitate to contact your local B. Braun Sales Representative for any questions regarding the B. Braun Product Portfolio or your local Kiro Grifols Sales Representative for questions related to KIRO® Oncology.

Yours sincerely,

B. Braun Melsungen AG

i.A. Loredana Tamburello



Senior Product Manager
Infusion Solutions – Ecoflac® plus
Phone: +49 151 40783859
e-mail: loredana.tamburello@bbraun.com

Kiro Grifols, S.L.

Naiara Telleria



Pharmacy Senior Manager
Phone: +34 943 252 349
e-mail: naiara.telleria@grifols.com

